

药物基因组学在精准医疗中的研究与应用



蒋正欣, PhD

先声诊断 高级医学经理

9.21.2019

精准医疗成败在于：精准用药

为何要做药物基因组学检测？

如何使患者精准用药获益最大化？

为何选择先声诊断和 MassARRAY® 平台？

精准医疗成败在于：精准用药

最新研究进展：临床证据与指南

重要临床获益性研究列举

先声诊断和 MassARRAY® 平台优势

为何药物反应存在显著个体化差异？

➤ 药代动力学 (PK) :

- 药物在体内吸收、分布、代谢和清除的动态规律

➤ 药效动力学 (PD) :

- 药物对机体的作用及其机制, 包括疗效和不良反应

➤ 药物基因组学 (PGx) 检测的目的:

- 明确引起个体药物反应差异的基因多态性
- 指导临床合理用药
- 疗效最大化 + 不良反应最小化
- 减少药物治疗费用

Drug	
Pharmacokinetic mechanisms	
CYP2B6	Efavirenz
CYP2C19	Clopidogrel, SSRIs, TCAs, voriconazole, proton pump inhibitors*
CYP2C9	Celecoxib*, phenytoin, warfarin
CYP2D6	Codeine, oxycodone, tramadol, SSRIs, TCAs, ondansetron, tamoxifen, atomoxetine
CYP3A5	Tacrolimus
DPYD	5-fluorouracil, capecitabine, tegafur
TPMT and NUDT15	Azathioprine, mercaptopurine, thioguanine
SLCO1B1	Simvastatin
UGT1A1	Atazanavir
Pharmacodynamic mechanisms	
CFTR	Ivacaftor
CYP4F2	Warfarin
G6PD	Rasburicase
HLA-B	Abacavir, allopurinol, carbamazepine, phenytoin
IFNL3 (IL28B)	Interferon
RYR1 and CACNA1S	Inhaled anesthetics, succinylcholine
VKORC1	Warfarin

SSRI=selective serotonin reuptake inhibitor. TCA=tricyclic antidepressant.
*Guidelines in progress.


Table 1: Drugs and genes with guidelines from the Clinical Pharmacogenetics Implementation Consortium for use in clinical practice

99 种药物 + 132 个 PGx 用药剂量指南 + 309 种药物的标签有 PGx 信息



根据 GRADE 和 AGREE II 指南制定标准, 推荐临床优先采用 **CPIC** 和 **CPNDS** 指南

Guideline	Scope and purpose (%)	Stakeholder involvement (%)	Rigor of development (%)	Clarity of presentation (%)	Applicability (%)	independence (%)	Overall score	Recommendation for use
CPIC	79.6	44.4	77.8	79.6	47.2	86.1	6/7	Recommended
RNPGX	66.7	29.6	43.1	72.2	33.3	38.9	3/7	Recommended with modifications
CPNDS	85.2	53.7	68.8	77.8	43.1	77.8	6/7	Recommended
DPWG	75.9	40.7	54.2	83.3	43.1	58.3	4/7	Recommended with modifications
Median (range)	77.8 (66.7–85.2)	42.6 (29.6–53.7)	61.5 (43.1–77.8)	78.7 (72.2–83.3)	43.1 (33.3–47.2)	68.1 (38.9–86.1)		




24种最重要的药物基因检测有指南或药物标签支持改变临床用药方案

VIP	Clinical Genetic Testing ^{1,2}	Drug (Guideline/Drug Label Organizations ³)	Clinical Impact	
<i>BRCA1</i>	Genetic testing for <i>BRCA1</i> mutations	Olaparib, rucaparib (FDA drug label)	Targeted treatment specific to genetic status	
<i>BRCA2</i>	Genetic testing for <i>BRCA2</i> mutations	Olaparib, rucaparib (FDA drug label)	Targeted treatment specific to genetic status	
<i>CACNA1S</i>	Genetic testing for <i>CACNA1S</i> mutations	Desflurane, enflurane, halothane, isoflurane, methoxyflurane, sevoflurane, succinylcholine (CPIC [2])	Alternate choice of medication to prevent serious ADR (risk of death)	
<i>CFTR</i>	Genetic testing for presence of <i>CFTR</i> G551D, F508del variants (+32 other variants now approved—found on ivacaftor drug label)	Ivacaftor (CPIC [7]), lumacaftor (when in formulation with ivacaftor) (FDA drug label)	Targeted treatment specific to genetic status	
<i>CYP2C19</i>	Genetic testing for presence of increased and decreased function alleles	Clopidogrel (DPWG [3,4], CPIC [10]) Amitriptyline, clomipramine, doxepin, imipramine, trimipramine (CPIC [32]—all tricyclic antidepressants listed, DPWG [3,4])—only imipramine Citalopram, escitalopram, sertraline (CPIC [31], DPWG [3,4]) Voriconazole (CPIC [6], DPWG [3,4])	Dosing adjustment/alternate choice of medication (risk of poor efficacy/ADRs)	
		Lansoprazole, omeprazole, pantoprazole (DPWG [3,4])	Increase attention/monitoring dose	
<i>CYP2C9</i>	Genetic testing for presence of decreased function alleles	Phenytoin (CPIC [20], DPWG [3,4])	Dosing adjustment to prevent serious ADR	
	Genetic testing for presence of decreased function alleles	Warfarin (CPIC [9], CPNDS [11], DPWG [3])	Dosing adjustment for optimal efficacy (avoiding excessive bleeding/clotting)	
<i>CYP2D6</i>	Genetic testing for presence of increased and decreased function alleles (recommendation may be based on genotype activity score)	Amitriptyline, also likely applicable to other TCAs: Clomipramine, desipramine, doxepin, imipramine, nortriptyline, trimipramine (CPIC [32]—as listed, DPWG [3,4]—only amitriptyline, clomipramine, doxepin, imipramine, nortriptyline) Aripiprazole, haloperidol, pimozide, zuclopenthixol (DPWG [3,4]), fluvoxamine (CPIC), paroxetine (CPIC [31]—both SSRIs listed, DPWG [3,4]—only paroxetine) Venlafaxine (DPWG [3,4]) Codeine (CPIC [23], DPWG [3,4], CPNDS [24]), tramadol (DPWG [3,4]) Flecainide, propafenone (DPWG [3,4]) Metoprolol (DPWG [3,4]) Tamoxifen (CPIC [26], DPWG [3,4], CPNDS [29]) Eliglustat (DPWG [3]) Tetrabenazine (FDA drug label)	Dosing adjustment/alternate choice of medication (risk of poor efficacy/ADRs)	
		Ondansetron, tropisetron (CPIC [8])	Alternate choice of medication to reduce risk of poor efficacy for UMs	
		Genetic testing for presence of two decreased function alleles	Atazanavir (CPIC [13])	Dosing adjustment to prevent serious ADR (jaundice)
		Genetic testing for presence of <i>UGT1A1</i> *1,*28,*36,*37 variants	Irinotecan (DPWG [3,4], French Group of Clinical Oncopharmacology (GPCO-Uncancer) & National Pharmacogenetics Network (RNPx) [30])	Dosing adjustment to prevent serious ADR (hematological/gastrointestinal toxicity)
		Genetic testing for presence of <i>UGT1A6</i> *4 (rs17863783) variant	Daunorubicin, doxorubicin (CPNDS [28])	Pediatric patients: Dosing adjustment to prevent serious ADR (cardiotoxicity)
		Genetic testing for homozygous <i>VKORC1</i> rs9934438 status	Acenocoumarol, phenprocoumon (DPWG [3,4])	Increase attention/monitoring dose
		Genetic testing for presence of <i>VKORC1</i> rs9923231 variant	Warfarin (CPIC [9], CPNDS [11], DPWG [3])	Dosing adjustment for optimal efficacy (avoiding excessive bleeding/clotting)

		Atomoxetine (CPIC [33], DPWG [3,4])	Increase attention/monitoring dose
<i>CYP3A5</i>	Genetic testing for presence of “normal” function and decreased function alleles	Tacrolimus (CPIC [16], DPWG [3,4])	Dosing adjustment to reduce risk of poor efficacy
		Warfarin (CPIC [9])	Dosing adjustment for optimal efficacy (avoiding excessive bleeding/clotting)
<i>CYP4F2</i>	Genetic testing for presence of <i>CYP4F2</i> *3 allele		
<i>DPYD</i>	Genetic testing for presence of decreased function alleles (recommendation based on genotype activity score)	Capecitabine, fluorouracil, tegafur (CPIC [25]—only capecitabine and fluorouracil, DPWG [3,4]—all three anti-neoplastics listed)	Dosing adjustment/alternate choice of medication (risk of ADR—death)
<i>DMD</i>	Genetic testing for presence of <i>DMD</i> mutation that is amenable to exon 51 skipping	Eteplirsen (FDA drug label)	Targeted treatment specific to genetic status
<i>F5</i>	Genetic testing for <i>F5</i> alleles	Estrogen-containing hormonal contraceptives (DPWG [3,4])	Alternate choice of contraceptive method to prevent serious ADR (venous thrombo-embolism)
<i>G6PD</i>	Genetic testing for presence of decreased function (class I, II, or III) alleles [x-linked—males 1 allele, females—2 alleles; if ambiguous result or female heterozygote—enzymatic testing to confirm activity levels]	Rasburicase (CPIC [18]) Pegloticase (FDA drug label, European Medicines Agency drug label) Primaquine (FDA drug label)	Alternate choice of medication to prevent serious ADR (acute hemolytic anemia)
		Carbamazepine (CPIC [21], CPNDS [22])	Dosing adjustment to prevent serious ADR (SCAR)
		Carbamazepine (CPIC [21], CPNDS [22]), phenytoin (CPIC [20]), oxcarbazepine (CPIC [21])	Dosing adjustment to prevent serious ADR (SCAR)
<i>HLA-A</i>	Genetic testing for presence of <i>HLA-A</i> *31:01 variant	Abacavir (CPIC [12], DPWG [3,4])	Dosing adjustment/alternate choice of medication (risk of poor efficacy/ADR—SCAR)
		Allopurinol (CPIC [17], American College of Rheumatology [19])	Dosing adjustment to prevent serious ADR (SCAR)
<i>HLA-B</i>	Genetic testing for presence of <i>HLA-B</i> *57:01 variant		
		Genetic testing for presence of <i>HLA-B</i> *58:01 variant	
<i>IFNL3</i>	Genetic testing for presence of <i>IFNL3</i> (IL28B) variant (rs12979860)	Peginterferon alfa-2a, peginterferon alfa-2b, ribavirin (CPIC [14])	Anticipated efficacy—consider in context of SDM and likely side effects
<i>POLG</i>	Mitochondrial genetic testing for <i>POLG</i> mutations	Divalproex sodium (FDA drug label, Health Canada/Santé Canada drug label)	Alternate choice of medication to prevent serious ADR (acute liver failure and death)
<i>RARG</i>	Genetic testing for presence of <i>RARG</i> rs2229774 variant	Daunorubicin, doxorubicin (CPNDS [28])	Pediatric patients: Dosing adjustment to prevent serious ADR (cardiotoxicity)
<i>RYR1</i>	Genetic testing for <i>RYR1</i> mutations	Desflurane, enflurane, halothane, isoflurane, methoxyflurane, sevoflurane, succinylcholine (CPIC [2])	Alternate choice of medication to prevent serious ADR (risk of death)
<i>SLCO1B1</i>	Genetic testing for presence of C allele at <i>SLCO1B1</i> rs4149056	Simvastatin (CPIC [5], DPWG [3])	Dosing adjustment/alternate choice of medication to prevent serious ADR (myopathy)
<i>TPMT</i>	Genetic testing for presence of decreased function alleles	Azathioprine, mercaptopurine, thioguanine (CPIC [15], DPWG [3,4])	Dosing adjustment/alternate choice of medication (risk of poor efficacy/ADRs)
		Genetic testing for presence of <i>TPMT</i> *2, *3A, *3B, *3C alleles	Cisplatin (CPNDS [27])

药物	基因	PGx 指南
→ 氯吡格雷	<i>CYP2C19</i>	CPIC, DPWG
华法林	<i>CYP2C9, VKORC1, CYP4F2</i>	CPIC, DPWG, CPNDS
辛伐他汀	<i>SLCO1B1</i>	CPIC, DPWG



CYP2C19与高危人群 (PCI-stent) 的主要心脏副作用相关, 与低危人群无显著相关性

Ref	Year	Ethnic	Population studied	n	PCI-stent (%)	Follow up	Endpoint	Polymorphisms	Outcomes (LOF vs no LOF)
High-risk patients (PCI-stent)									
Collet (29)	2009	Europeans	ACS	259	86	>4 years	MACE (CV death, ACS, urgent PCI)	CYP2C19*2	HR 5.38 (2.32-12.47) p ≤ 0.0001
							ST definite	CYP2C19*2	HR 6.04 (1.75-20.80) p = 0.004
Mega (22)	2009	84% Europeans	ACS stent (TRITON)	1477	100	15 months	MACE (CV death, ACS, stroke)	CYP2C19*2	HR 1.53 (1.07-2.19) p = 0.01
							ST definite	CYP2C19*2	HR 3.09 (1.19-8.0) p = 0.02
Mega (17)	2010	84% Europeans	ACS stent (TRITON)	2905	100	15 months	MACE (CV death, ACS, stroke)	CYP2C19*2 and ABCB1	ABCB1 TT vs CT/CC: HR 1.72 (1.22-2.44) p = 0.002 CYP2C19*2 + ABCB1 HR 1.97 (1.38-2.82) p = 0.0002
Simon (30)	2009	Europeans	ACS	2208	68.7	12 months	MACE (death any cause, ACS, stroke)	CYP2C19 and ABCB1	CYP2C19: HR 1.98 (1.10-3.58) ABCB1: HR 1.72 (1.20-2.47)
Sorich (33)	2010	84% Europeans	ACS stent (TRITON)	13608	100	15 months	MACE (CV death, ACS, stroke)	CYP2C19 LOF	OR 1.63 (1.45-1.81) p < 0.0001
Shuldiner (31)	2009	Europeans	PCI	227	100	12 months	MACE (CV death, ACS, stroke, PCI)	CYP2C19*2	HR 2.42 (1.18-4.99) p = 0.02
Wallentin (35)	2010	Europeans	ACS	10285	60	12 months	MACE (CV death, ACS, stroke)	CYP2C19	HR at 30 days: p = 0.028
								CYP2C19 and ABCB1	HR 1.2 (1.0-1.4) p = 0.047**
Low-risk patients									
Pare (53)	2010	Europeans-latin american	ACS stable	5059	14.5	12 months	MACE (CV death, ACS, stroke)	CYP2C19*2	p = 0.32

Data are showed as: OR: odds ratio, HR: hazard ratio, (95%CI), p-value

CYP2C19与高危人群的主要心脏副作用相关, CPIC和DPWG指南建议高危人群检测

Ref	Year	Ethnic	Population studied	n	PCI-stent (%)	Follow up	Endpoint	Polymorphisms	Outcomes (intervention group vs control group)
Non RCT									
Sánchez-Ramos (37)	2016	Europeans	ACS-PCI-stent	719	100	1 year	MACE (CV death, ACS, stroke)	CYP2C19*2, *3 and ABCB1	HR 0.63 (0.41-0.97) p=0.037
							ST definite	CYP2C19 *2, *3 and ABCB1	HR 1.27 (0.08-20.2) p=0.87
							Urgent revascularization*	CYP2C19 *2, *3 and ABCB1	HR 0.63 (0.31-1.28) p=0.20
RCT									
→ Shen (39)	2016	Asians	CAD-PCI	628	100	1 month 6 months 12 months	MACE (composite of death from any cause, myocardial infarction, or target vessel revascularization)	CYP2C19*2	1.3% vs. 5.6%, P = 0.003 3.2% vs. 7.8%, P = 0.012 4.2% vs. 9.4%, P = 0.010
Roberts (Rapid Gene) (40)	2012	Europeans	ACS or stable angina/ stent	187	100	7 days	high on-treatment platelet reactivity	CYP2C19*2	0% vs 30% p=0.0092
Roberts (RAPID STEMI study) (41)	2016	Europeans	STEMI- stent	102	100	1 month	high on-treatment platelet reactivity	CYP2C19*2, *17 and ABCB1 TT	OR=0.15 p=0.03
→ Xie (42)	2013	Asians	CAD-PCI	600	100	180 days	MACE (death from any cause, MI, stroke, ischemia)	CYP2C19*2,*3	1.0% and 6.2%, P<0.01
Notarangelo Pharmclo (43)	2018	Europeans	ACS	888	No data	12 months	MACE (CV death, nonfatal IM, nonfatal stroke)	CYP2C19*2, *17 and ABCB1	HR 0.58 (0.43-0.78) p<0.001
Bergmeijer (Popular genetics) (44)	ongoing	Europeans	STEMI-stent	2500	100	15 months	MACE (CV death, ACS, stroke)	CYP2C19*2, *3	
Tailor-PCI (NCT01742117)	ongoing	Europeans	ACS or CAD/ stent	5000	100	12 months	MACE (non-fatal MI, non-fatal stroke, severe recurrent ischemia, CV death, and ST)	CYP2C19*2,*17	

Data are showed as: OR: odds ratio, HR: hazard ratio, (95%CI), p-value

CPIC 和 DPWG 指南建议检测 CYP2C19 指导氯吡格雷用药方案

PGx指南	基因代谢型	风险描述	适用人群	治疗建议
CPIC	CYP2C19 EM	血小板抑制和聚集正常	ACS-PCI	按照药物标签常规使用
	CYP2C19 IM	血小板抑制减弱、血小板聚集和心血管不良反应风险升高	ACS-PCI	选择替代药物 (例如: 普拉格雷、替格瑞洛)
	CYP2C19 PM	血小板抑制极大减弱、血小板聚集和心血管不良反应风险升高	ACS-PCI	选择替代药物 (例如: 普拉格雷、替格瑞洛)
	CYP2C19 UM	血小板抑制增强、血小板聚集减弱	ACS-PCI	按照药物标签常规使用
DPWG	CYP2C19 IM	活性代谢物减少, 严重心脑血管不良反应风险升高	PCI, Stroke, TIA	选择替代药物 (例如: 普拉格雷、替格瑞洛) 或增加剂量至150 mg/d (初始剂量600 mg)
	CYP2C19 PM	活性代谢物减少, 严重心脑血管不良反应风险升高	PCI, Stroke, TIA	选择替代药物 (例如: 普拉格雷、替格瑞洛)
	CYP2C19 UM	活性代谢物增加, 严重心脑血管不良反应风险可能降低, 出血风险可能升高	普适	按照药物标签常规使用

CPIC 2013 update: Scott SA, et al., Clin Pharmacol Ther. 2013

DPWG 2018 update

CYP2C19, VKORC1, CYP4F2 基因检测指导华法林剂量调整:


- 可改善治疗范围时间比例 (%TTR)
- 降低出血和静脉血栓风险

Ref	Year	Ethnic	Population studied	n	Follow up	Endpoint	Polymorphisms	Homozygotes action required	Outcomes (intervention group vs control group)	Availability test	Dose in non-genotype group
Pirmohamed (EU-PACT) (3)	2015	2% non-European	AF (72.1%) VT (27.9%)	455	12 weeks	%TTR	CYP2C9*2,*3 VKORC1	VKORC1: 17% CYP2C9*2 and *3: 3.4%	67.4% vs 60.3%, p < 0.001	2h	Fixed-dose strategy
Kimel (COAG) (79)	2015	33% non-European	AF (23%) DVT or PE (56%)	1015	4 weeks	%TTR	CYP2C9*2,*3 VKORC1	VKORC1: 11% CYP2C9*2 and *3: 1%	45.2% vs 45% p = 0.91	Not before the 1 st dose for 55% of patients	Clinical dosing algorithm
Gage (GIFT) (82)	2017	91% European	Hip or knee arthroplasty	1650	30 and 60 days	Composite (major bleeding, INR≥4, VT, death)	CYP2C9*2,*3 CYP4F2	NA	RR 0.73 (0.56-0.95) p = 0.02	NA	NA


AF: atrial fibrillation, DVT: Deep-vein thrombosis, PE: pulmonary embolism. %TTR: percentage time in therapeutic range. NA: Not applicable

CPIC, DPWG, CPNDS 指南建议基因检测指导华法林剂量调整

成人



儿童



Warfarindosing.org


Gage, IWPC 算法

Johnson JA, et al., Clin Pharmacol Ther. 2017

Gage BF, et al., Clin Pharmacol Ther. 2008



International Warfarin Pharmacogenetics Consortium, N Engl J Med. 2009

药物	基因	PGx 指南
氯吡格雷	<i>CYP2C19</i>	CPIC, DPWG
华法林	<i>CYP2C9, VKORC1, CYP4F2</i>	CPIC, DPWG, CPNDS
→ 辛伐他汀	<i>SLCO1B1</i>	CPIC, DPWG



SLCO1B1*5 与肌肉病变不良反应风险相关

STRENGTH study (n=452: 99 cases+353 controls)




Brunham study (n=108: 25 cases+83 controls)

Table 2 Association between rs4149056 genotype and myopathy

Statin-type	Cases			Controls			C allele frequency		
	T/T	T/C	C/C	T/T	T/C	C/C	Cases	Controls	P value
All	15	8	2	57	20	6	0.24	0.19	0.21
Simva	5	6	1	27	10	2	0.33	0.18	0.042
Atorva	7	2	1	24	7	3	0.20	0.19	0.48

SLCO1B1*5 (rs4149056: C allele)



Voora D, et al., J Am Coll Cardiol. 2009

Brunham LR, et al., Pharmacogenomics J. 2012

CPIC 和 DPWG 指南建议检测 *SLCO1B1* 指导辛伐他汀用药方案

PGx指南	基因多态性	风险描述	治疗建议
CPIC	rs4149056: TT	肌肉病变风险较低	按适应症指南选择药物剂量常规使用
	rs4149056: CT	肌肉病变风险适中	降低剂量或选择替代药物（例如：普伐他汀、瑞舒伐他汀），并常规监测肌酸激酶
	rs4149056: CC	肌肉病变风险较高	降低剂量或选择替代药物（例如：普伐他汀、瑞舒伐他汀），并常规监测肌酸激酶
DPWG	rs4149056: CT	肝细胞摄取减少，血药浓度和肌肉病变风险升高	<ol style="list-style-type: none"> 1. 选择替代药物（例如：普伐他汀、瑞舒伐他汀、氟伐他汀） 2. 若无替代药物则最大剂量不超过40 mg/d，并在出现肌肉症状时及时就医
	rs4149056: CC	肝细胞摄取减少，血药浓度和肌肉病变风险升高	选择替代药物（例如：普伐他汀、瑞舒伐他汀、氟伐他汀），并考虑他汀诱导肌肉病变的其他风险因子

Mayo Clinic	时效性检测 (Reactive)	预判性检测 (Preemptive)
检测时机	用药时	用药前 / 出生时
检测目的	指导当前特定用药 / 剂量	指导未来用药相关处方
检测基因	1 个或多个基因	多基因 Panel / 全基因组
结果参考	可能不再参考	多次参考
检测时效	要求高 / 越快越好	要求不高 / 随用随取
检测项目	华法林 CYP2C9 / VKORC1 氯吡格雷 CYP2C19 他汀类降血脂药 SLCO1B1 / APOE	高血压用药 Panel 精神用药 Panel 冠心病联合用药 Panel 心血管用药 Panel CPIC 指南用药 Panel 综合型 FDA 标签用药 Panel 全基因组用药 Panel

先声诊断：用精准医疗帮助每个人

Live longer, live better



- 先声诊断是中国领先的医药集团先声控股旗下**精准医疗解决方案提供商**
- 致力于在人类健康的重大挑战领域创造**精准、全面、快捷、可及**的诊断产品和服务体系

临床解决方案



面向医务工作者和医疗机构


通过整合多种领先的分子检测平台和生物信息学分析能力，在肿瘤、中枢神经系统、感染、药物基因组学、自身免疫疾病等治疗领域提供精准医疗的决策基础

个人检验服务



面向患者和健康人群

通过深度整合科技、数据及制造能力提供基于自采样、寄样的专业医学检验服务以创新的服务模式提升公共医疗的质量、效率和体验



- 先声诊断与美国Agena达成战略合作，共同开拓核酸质谱平台MassARRAY在中国伴随诊断和药物基因组学领域的产业化应用
- **独家权益：** MassARRAY机器 + 药物基因组学PGX试剂盒的NMPA注册及销售
- **其它权益：** 肿瘤、液体活检、甲基化及病原体耐药试剂盒的NMPA注册及销售


- 先声诊断与罗氏诊断达成战略合作，共同致力于**新一代测序产品**的开发与应用，为更多终端用户提供高品质的检测结果，满足转化医学精准治疗的细致需求
- 借助罗氏诊断产品开发能力，**完善先声的分子诊断产品线**，为临床诊断提供更加符合未来发展的解决方案

- 先声诊断与安捷伦达成战略合作
- 双方将发挥各自的优势，紧密合作，共同推动精准医疗发展，帮助科学家和医生应对人类健康的重大挑战



肿瘤

- 覆盖**实体瘤**和**血液肿瘤**伴随诊断服务 - 患者全生命周期管理式诊断
- 遗传/易感检测
- 早期筛查
- 分子分型/用药指导
- 预后评估
- 实时动态监控



中枢神经

- 针对**脑部肿瘤**诊疗产品, 对不同基因型的患者提供差异化治疗方案
- 针对6大**神经系统自身免疫性疾病**, 进行自身抗体的全面检测, 提供疑难疾病的综合诊断



感染病原体

- 提供**感染病原体**快速分析和诊断的全套解决方案
- 用药/耐药/药敏监测
- 急重症/复杂感染检测
- 普检
- 初筛
- 家庭自检



药物基因组

- 指导临床选择**最佳治疗药物**和**最适剂量**提高药物疗效、避免不良反应、改善预后、节约医疗成本, 实现个体化治疗
- 对FDA批准、CPIC推荐并在国内上市的药物基本全覆盖

先声诊断 - 国内该技术平台在药物基因组领域产品类型最丰富的公司



药物不良反应频发

- 2017年全国药品不良反应事件报告约**143万份**
- 药物不良反应已成为我国居民死亡的**第四大原因**



遗传是个体化差异首要原因

- 约**60%**的药物反应个体化差异是遗传因素决定的



用药检测是重要解决方法

- ✓ 评估药物疗效和毒副作用，选择合适药物
- ✓ 调整用药剂量，减少过度治疗或治疗不足
- ✓ 避免药物不良反应发生

基于MassARRAY®平台的药物基因组检测组合

心脑血管疾病
个体化用药检测

降脂药物
个体化用药检测

降压药物
个体化用药检测

抗血栓药物
个体化用药检测

神经与精神疾病
个体化用药检测

儿童
安全用药检测

产品 | 神经与精神疾病个体化用药检测



全面: 5 大类神经与精神疾病 + 52 种常见药物 + 23 个基因 + 30 个多态性位点

专业: 基于国内外最新权威指南和数据库

细致: 3 个小套餐 - 抗癫痫药、抗抑郁与焦虑药、抗精神病药

	药物标签				PGx指南			
	FDA	HCSC	PMDA	EMA	CPIC	DPWG	CPNDS	中国卫计委
卡马西平	HLA-B, HLA-A	HLA-B, HLA-A	HLA-B, HLA-A		HLA-B, HLA-A		HLA-B, HLA-A	HLA-B
奥卡西平	HLA-B				HLA-B			
苯妥英	HLA-B, CYP2C9	HLA-B			HLA-B, CYP2C9	CYP2C9		HLA-B, CYP2C9
氟巴占	CYP2C19							
丙戊酸	POLG	POLG						
地西洋	CYP2C19							
西酞普兰	CYP2C19	CYP2C19			CYP2C19	CYP2C19		
艾司西酞普兰	CYP2C19		CYP2C19		CYP2C19	CYP2C19		
氟西汀	CYP2D6							
氟伏沙明	CYP2D6				CYP2D6			
丙咪嗪	CYP2D6				CYP2C19, CYP2D6	CYP2C19, CYP2D6		CYP2D6
氯丙咪嗪	CYP2D6				CYP2C19, CYP2D6	CYP2D6		CYP2D6
舍曲林					CYP2C19	CYP2C19		
阿米替林	CYP2D6				CYP2C19, CYP2D6	CYP2D6		CYP2C19, CYP2D6
地昔帕明	CYP2D6				CYP2D6			CYP2D6
去甲替林	CYP2D6	CYP2D6			CYP2D6	CYP2D6		CYP2D6
帕罗西汀	CYP2D6				CYP2D6	CYP2D6		
文拉法辛	CYP2D6					CYP2D6		
多塞平	CYP2C19, CYP2D6				CYP2C19, CYP2D6	CYP2D6		
普罗替林	CYP2D6							
曲米帕明	CYP2D6				CYP2C19, CYP2D6			
沃替西汀	CYP2D6	CYP2D6		CYP2D6				
度洛西汀	CYP2D6			CYP2D6				
氟哌啶醇						CYP2D6		
阿立哌唑	CYP2D6	CYP2D6		CYP2D6		CYP2D6		
氯氮平	CYP2D6							
奥氮平	CYP2D6			CYP2D6				
利培酮	CYP2D6	CYP2D6						
奋乃静	CYP2D6		CYP2D6					
珠氯噻醇						CYP2D6		
硫利达嗪	CYP2D6							
阿托莫西汀	CYP2D6	CYP2D6	CYP2D6		CYP2D6	CYP2D6		

全面: 4 大类风湿免疫疾病 + 19 种常见药物 + 29 个基因 + 41 个多态性位点

专业: 基于国内外最新权威指南和数据库

细致: 3 个小套餐 - 改善病情的抗风湿药、生物制剂和抗痛风药

	药物标签			PGx指南			
	FDA	HCSC	PMDA	CPIC	DPWG	ACR	中国卫计委
硫唑嘌呤	NUDT15, TPMT	TPMT	TPMT	NUDT15, TPMT	TPMT		TPMT
他克莫司				CYP3A5	CYP3A5		CYP3A5
羟氯喹	G6PD						
柳氮磺吡啶	G6PD, NAT2	G6PD	G6PD				
别嘌醇			HLA-B	HLA-B		HLA-B	HLA-B

全面: 9大疾病领域 + 140种儿童常用药物 + 33个基因 + 50个多态性位点

专业: 基于国内外最新权威指南和数据库

细致: 2种端口 - 面向临床、争对患者

适用人群	人群特征	检测意义
患病儿童	发热疼痛患者 哮喘咳嗽患者 神经与精神病患者 心脑血管疾病患者 感染病患者 糖尿病患者 消化系统疾病患者 麻醉病人 免疫化疗患者	精准用药: 从剂量、药效、不良反应等多方面给予个性化用药提示, 改善患者疗效并减少或避免严重副作用
健康儿童	关注用药风险儿童 曾有过药物不良反应史的儿童 家族成员中有过药物不良反应史的儿童	风险评估: 评估儿童用药风险 指导安全用药: 根据风险评估结果, 选择合适药物, 有效规避潜在用药不当风险

MassARRAY临床应用




药物基因组全线产品



脑胶质瘤精准检测



泛癌种化疗用药



遗传性耳聋检测



营养元素基因检测



叶酸代谢基因检测

更多出生缺陷及遗传领域临床应用开发中

先声诊断与Agena达成战略合作

Agena授权了合作开发及NMPA注册权益





MassARRAY 系统将质谱分析、高灵敏度和稳定性的化学试剂与先进的数据分析软件相结合
以满足基因组学实验室的检测设计、验证和性能需求

高灵活性

简化了在实验室中设计和实施开展新测试项目的流程
PCR反应条件均一化
检测试剂通用化
位点探针设计便捷

高性价比

检测成本有明显的价格优势
无需荧光标记、低成本定量分析

高精度度

MassARRAY[®]核酸质谱仪检测技术灵敏度高，检测结果准确
FDA将其特许作为第二代测序仪
illumina MiSeqDx 的验证平台

多功能性


可检测SNP/Indel基因分型、
基因融合、染色体结构变异、
DNA甲基化、CNV等各类已知
变异

高通量性

在单个反应孔中实现
高度多重性的能力（高达40重）
可通宵运行，最大程度提高通量

全自动化

从检测到结果报告全程8小时
无需生信分析
点样及质谱飞行过程全自动化
最大程度降低污染风险




打包第三方检验服务

按相应政策法规为医院提供符合质量规范的临检服务及医学支持



技术供应商

基础项目本地化/人才培养/PCR实验室建设等



区域检验中心（暨国重分中心）建设

国家重点实验室联手打造集约化区域检验中心，提升临检与科研发展水平

Thank You



别嘌醇、卡马西平、华法林和氯吡格雷在临床进行基因检测经济有效

	参考研究	适应症	标记物人群频率	有效性	经济性
阿巴卡韦 HLA-B*57:01 ADR: 过敏	Hughes et al. 2004	HIV	UK: 3.7%	Test: 1.9% ADR; No test: 3.7% ADR	\$490/ADR avoided for testing
	Schackman et al. 2008	HIV	USA: 5.7%	Test: 0% ADR, No test: 2.73% ADR	Test dominates other strategies
	Kapoor et al. 2014	HIV	Han Chinese 1.1%	Test vs No test: Improve quality of life	Not cost-effective for testing
	Calatrava et al. 2010	HIV	Spain: 5.6%	Test: 4.25% ADR; No test: 7.8% ADR	\$630/ADR avoided for testing
别嘌醇 HLA-B*58:01 ADR: 过敏	Park et al. 2015	Gout	Korea: 12.2%	Test: 0% ADR, No test: 2.19% ADR	Test dominates other strategies
卡马西平 HLA-B*15:02, HLA-A*31:01 ADR: SJS, TEN	Dong et al. 2012	Epilepsy	Singapore: 14.87%	Test vs No test: Improve quality of life	Test dominates valproate
华法林 CYP2C9*2, *3, VKORC1 ADR: 出血	Leey et al. 2009	AF	USA: NA	Test vs No test: reduce major bleeding	Test dominates other strategies
	You et al. 2009	Anticoagulation	USA: NA	Test vs No test: reduce ADR	\$170192/ADR averted for testing
	Schalekamp et al. 2006	Thromboembolic	Netherlands: 36%	Test: 4.6% ADR, No test: 5.1% ADR	\$4233/ADR averted for testing
	You et al. 2012	AF	USA: NA	Test vs No test: Improve quality of life	Test dominates standard care
	You et al. 2004	Thromboembolic	USA: 36%	Test: 9.58% ADR, No test: 10.48% ADR	\$5778/ADR averted for testing
氯吡格雷 CYP2C19 ADR: 主要心血管不良反应	Lala et al. 2012	ACS	USA: 27.1%	Test vs No test: Improve quality of life	Test dominates other strategies
	Panattoni et al. 2012	ACS	Asian 29%	Test vs No test: Improve quality of life	Test dominates prasugrel
	Reese et al. 2012	ACS, MI, stroke	USA: 27%	Test vs No test: Improve quality of life	Test dominates other strategies
	Sorich et al. 2013	ACS	Australia: 28.2%	Test vs No test: Improve quality of life	Test dominates ticagrelor