
# Phenotype MicroArrays

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


**Results**  
**Behavior on Biolog Phenotype MicroArray Plates.** We tested our parental strain NCM3722 and strains with mini Tn5 insertions in several genes of the b1012 operon on Biolog (Hayward, CA) Phenotype MicroArray plates (PM3) for nitrogen catabolism, which couple utilization of various N sources to reduction of a tetrazolium dye and production of purple color (10, 11). We also tested a strain carrying an *ntrB*(Con) mutation, NCM3876, and corresponding mutant derivatives. Strains NCM3722 and NCM3876 [*ntrB*(Con)] could respire with uracil or uridine as N source at room temperature, whereas mutant derivatives could not (Fig. 2 and data not shown). Parental strain NCM3876 [*ntrB*(Con)] could also respire with thymidine or thymine as N source (data not shown). None of the strains could respire with pyrimidines as N source at 37°C.


Low, Kustu, and coworkers., PNAS 103, 5114(2006). <https://doi.org/10.1073/pnas.0600521103>

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
**Figure 2 | Evolution of bacterial metabolic growth phenotypes.** Genetic distances in the figure are based on bacterial 16S rRNA sequences. **a.** The evolution of phenotypic similarity in the usage of carbon sources for biomass synthesis. The colours represent the point density at a given genetic distance for all pairwise comparisons between metabolic models ( $n = 26,106$ ). The black line shows a three-parameter exponential fit to the computational predictions; the red line shows a moving average of the predictions. **b.** Like **a**, but phenotypic similarity in the usage of nitrogen sources for biomass production across metabolic models ( $n = 36,856$ ). **c.** Experimental analysis of the long-term phenotypic divergence trend. Grey points represent pairwise comparisons of

carbon usage phenotypes (Biolog data) between 40 bacterial species ( $n = 780$ ). The black squares represent the average values of experimental phenotypic similarity at different divergence distances. The black line represents an exponential fit to the experimental phenotypic similarity data; the red line represents an exponential fit to the computationally predicted phenotypic similarity data for biomass synthesis. **d.** The variability of experimental and computationally predicted phenotypic similarity at different divergence distances. The variability was quantified by the coefficient of variation, defined as the ratio of the standard deviation to the mean. Error bars in **c** and **d** represent s.e.m. obtained on the basis of 10,000 bootstrap re-samplings of the considered species.

**Figure 3 | Distribution of phenotypic similarity at different levels of bacterial taxonomic classification.** The distributions of phenotypic similarity in the usage of carbon sources for biomass synthesis were obtained based on computational simulations of metabolic models ( $n = 26,106$ ). The dashed line connects the average values (black dots) of computational predictions at each taxonomic level. The red asterisks in the figure indicate the average values of experimental data obtained using Biolog arrays (see Methods).


**Extended Data Figure 4 | Evolutionary relationship between 40 species for which experimental phenotypic microarray data were considered.** The cladogram is based on 16S bacterial rRNA sequences. Different colours indicate different bacterial clades. The tree was rooted in the *M. luteus* rRNA sequence.

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
Yuta Otsuka et al., Nucleic Acids Research 43, D606–D617(2015).

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
- 
- Figure 1: Metabolic pathways and carbon source utilization of *H. pylori*.**
- Carbon Source Utilization (AUC) Data:**
- | Carbon Source          | AUC (approx.) |
|------------------------|---------------|
| D-Ribose               | 9200          |
| L-Lyxose               | 8800          |
| L-Arabinose            | 8500          |
| D-Xylose               | 8200          |
| D-Arabinose            | 7800          |
| D-Mannose              | 6500          |
| D-Malic Acid           | 5800          |
| Pyruvic Acid           | 4800          |
| 2-Deoxy-D-Arabinose    | 4200          |
| Fumaric Acid           | 4100          |
| L-Asparagine           | 4000          |
| D-Alanine              | 4000          |
| L-Alanine              | 3900          |
| Succinic Acid          | 3800          |
| Bromosuccinic Acid     | 3500          |
| D-Gluconic Acid        | 3200          |
| D-Glucuronic Acid      | 2800          |
| D-Fructose-6-phosphate | 2500          |
| L-Lysine               | 2000          |
| D-Lysine               | 1500          |
| D-Galactose            | 1000          |
- Metabolic Pathways:**
- Glycolysis:**  $\alpha$ -D-glucose-1P  $\leftrightarrow$   $\alpha$ -D-glucose  $\leftrightarrow$   $\alpha$ -D-glucose-6P  $\leftrightarrow$   $\beta$ -D-glucose  $\leftrightarrow$   $\beta$ -D-glucose-6P  $\leftrightarrow$   $\beta$ -D-fructose-6P  $\leftrightarrow$   $\beta$ -D-fructose-1,6P<sub>2</sub>  $\leftrightarrow$  glyceraldehyde-3P  $\leftrightarrow$  dihydroxyacetone  $\leftrightarrow$  glycerate-1,3P<sub>2</sub>  $\leftrightarrow$  glycerate-3P  $\leftrightarrow$  glycerate-2P  $\leftrightarrow$  phosphoenolpyruvate  $\leftrightarrow$  pyruvate  $\leftrightarrow$  acetyl-coA.
  - Entner-Dooudoroff Pathway:** D-glucose-6P  $\leftrightarrow$  D-glucono-1,5-lactone-6P  $\leftrightarrow$  6-phospho-D-gluconate  $\leftrightarrow$  2-dehydro-3-deoxy-D-gluconate-6P  $\leftrightarrow$  pyruvate.
  - Non-oxidative Pentose Phosphate Pathway:**  $\beta$ -D-fructose-6P  $\leftrightarrow$  D-erythrose-4P  $\leftrightarrow$  D-sedoheptulose-7P  $\leftrightarrow$  D-ribose-5P  $\leftrightarrow$  D-xylulose-5P  $\leftrightarrow$  D-ribose.
  - TCA Cycle:** acetyl-coA  $\leftrightarrow$  citrate  $\leftrightarrow$  oxaloacetate  $\leftrightarrow$  L-glutamate  $\leftrightarrow$  isocitrate  $\leftrightarrow$  oxalosuccinate  $\leftrightarrow$   $\alpha$ -ketoglutaric acid  $\leftrightarrow$  succinyl-coA  $\leftrightarrow$  succinate  $\leftrightarrow$  fumarate  $\leftrightarrow$  malate  $\leftrightarrow$  oxaloacetate.
- Legend:**
- Core carbon source used by all *H. pylori* strains, as indicated by PM data
  - $\cdots \rightarrow$  Involves series of reactions
  - $\rightarrow$  or  $\leftrightarrow$  Enzymatic reactions that able to be performed by *H. pylori*

**Measured metabolic activity (Phenotype microarray)**

min max

**Predicted metabolic activity (Ensemble FBA)**

min max

**Inorganic compounds**

**Protein amino acids**

**Nitrogen bases**

Root19 - Measured  
Root6201 - Measured  
Root497 - Measured  
Root19 - Predicted  
Root6201 - Predicted  
Root497 - Predicted

Ammonia  
Hydroxylamine  
Nitrate  
Nitrite  
Glycine  
L-Alanine  
L-Arginine  
L-Asparagine  
L-Aspartic  
L-Cysteine  
L-Glutamic  
L-Glutamine  
L-Histidine  
L-Isoleucine  
L-Leucine  
L-Lysine  
L-Methionine  
L-Phenylalanine  
L-Proline  
L-Serine  
L-Threonine  
L-Tryptophan  
L-Tyrosine  
L-Valine  
Adenine  
Cytosine  
Guanine  
Thymine  
Uracil  
Xanthine

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



Utilization of different compounds as sole carbon source by *Halomonas* TD01

| Category         | Plate position | Carbon source         | Growth observed | Omnilog unit (OU) |
|------------------|----------------|-----------------------|-----------------|-------------------|
| Negative control | PM1-A01        | Without carbon source | -               | 68                |
| Carbohydrate     | PM1-A10        | D-Trehalose           | +               | 276               |
|                  | PM1-B02        | D-Sorbitol            | +               | 265               |
|                  | PM1-B03        | Glycerol              | +               | 257               |
|                  | PM1-B11        | D-Mannitol            | +               | 269               |
|                  | PM1-C07        | D-Fructose            | +               | 273               |
|                  | PM1-C09        | α-D-Glucose           | +               | 268               |
|                  | PM1-C10        | Maltose               | +               | 268               |
|                  | PM1-D09        | α-D-Lactose           | +               | 204               |
|                  | PM1-D11        | Sucrose               | +               | 260               |
|                  | PM1-A02        | L-Arabinose           | +               | 165               |
|                  | PM1-B08        | D-Xylose              | +               | 159               |
| Carboxylic acid  | PM1-F11        | D-Cellobiose          | +               | 189               |
|                  | PM1-A05        | Succinic acid         | +               | 269               |
|                  | PM1-B06        | D-Gluconic acid       | +               | 267               |
|                  | PM1-B09        | L-Lactic acid         | +               | 256               |
|                  | PM1-B10        | Formic acid           | +               | 138               |
|                  | PM1-C03        | DL-Malic acid         | +               | 265               |
|                  | PM1-C08        | Acetic acid           | +               | 258               |
|                  | PM1-F02        | Citric acid           | +               | 268               |
|                  | PM1-H08        | Pyruvic acid          | +               | 276               |


J. Zhang, et al., Journal of Biotechnology 316, 1-5 (2020)  
<https://doi.org/10.1016/j.jbiotec.2020.04.007>

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

PM plot of the utilization of Aga by E. coli O157:H7 EDL933 and EC4045, a spinach isolate. Data derived from the increase in intensity of the color of the reduced dye A, expressed in AU (maximum of 500 AU), are plotted against time (48 h) using the PM software from Biolog. E. coli O157:H7 EDL933 utilizes Aga as carbon and nitrogen sources, whereas EC4045, a spinach isolate, cannot.

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




- A**

2-Aminoethanol  
2-Deoxy-D-Ribose  
2'-Deoxyadenosine  
4-Hydroxybenzoic Acid  
Acetamide  
Acetic Acid  
Acetoacetic Acid  
Adenosine  
Butyric Acid  
Capric Acid  
Citric Acid  
D,L-Citramalic Acid  
D,L-Malic Acid  
D,L-a-Glycerol Phosphate  
D-Fructose  
D-Fructose-6-Phosphate  
D-Galactose  
D-Gluconic Acid  
D-Glucose-1-Phosphate  
D-Glucose-6-Phosphate  
D-Mannose  
D-Melibiose  
D-Raffinose  
D-Ribose  
D-Saccharic Acid  
D-Serine  
D-Sorbitol  
D-Trehalose  
Dihydroxyacetone  
Formic Acid  
Fumaric Acid  
Glycerol  
Glycine  
Glycolic Acid  
Glyoxylic Acid  
Hydroxy-L-Proline  
Inosine  
L-Alanine  
L-Arabinose  
L-Arginine  
L-Asparagine  
L-Aspartic Acid  
L-Glutamic Acid  
L-Glutamine  
L-Histidine  
L-Homoserine  
L-Isoleucine  
L-Lactic Acid  
L-Leucine  
L-Lysine  
L-Malic Acid  
L-Methionine  
L-Ornithine  
L-Phenylalanine  
L-Proline  
L-Serine  
L-Threonine  
L-Valine  
N-Acetyl-D-Glucosamine  
N-Acetyl-L-Glutamic Acid  
Oxalic Acid  
Propionic Acid  
Putrescine  
Pyruvic Acid  
Succinic Acid  
Sucrose  
Thymidine  
Uridine  
a-D-Glucose  
a-D-Lactose  
a-Hydroxybutyric Acid  
b-Hydroxybutyric Acid  
b-Hydroxypropyruvic Acid  
m-Inositol

**B**

Acetamide  
Adenine  
Adenosine  
Agmatine  
Allantoin  
Ammonia  
Cytidine  
Cytosine  
D,L-a-Amino-N-Butyric Acid  
D-Glucosamine  
D-Serine  
Ethanolamine  
Ethylamine  
Glycine  
Guanine  
Guanosine  
Histamine  
Inosine  
L-Alanine  
L-Arginine  
L-Asparagine  
L-Aspartic Acid  
L-Citrulline  
L-Cysteine  
L-Glutamic Acid  
L-Glutamine  
L-Histidine  
L-Homoserine  
L-Isoleucine  
L-Leucine  
L-Lysine  
L-Methionine  
L-Ornithine  
L-Phenylalanine  
L-Proline  
L-Serine  
L-Threonine  
L-Tryptophan  
L-Tyrosine  
L-Valine  
Methylamine  
N-Acetyl-D,L-Glutamic Acid  
N-Acetyl-D-Glucosamine  
Nitrate  
Nitrite  
Putrescine  
Thymidine  
Uracil  
Urea  
Uric Acid  
Uridine  
Xanthine  
Xanthosine

**C**

| Carbon sources          | TP | TN | FP | FN | Accuracy | Sensitivity | Specificity | Positive Predicted Value |
|-------------------------|----|----|----|----|----------|-------------|-------------|--------------------------|
| <b>Carbon sources</b>   |    |    |    |    |          |             |             |                          |
| Heterotrophy            | 38 | 18 | 8  | 10 | 0.76     | 0.79        | 0.69        | 0.83                     |
| <b>Nitrogen sources</b> |    |    |    |    |          |             |             |                          |
| Heterotrophy            | 27 | 17 | 8  | 1  | 0.83     | 0.96        | 0.68        | 0.77                     |

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